

Clinical perspective

Measurement of toe systolic pressures: a technique paper

Tehan P, Fox M and Mills JL

Keywords vascular assessment, toe systolic pressure, toe-brachial index, technique, podiatry, woundcare**For referencing** Tehan P et al. Measurement of toe systolic pressures: a technique paper. *Wound Practice and Research* 2021; 29(3):148-153.**DOI** <https://doi.org/10.33235/wpr.29.3.148-153>

Submitted 2 June 2021, Accepted 30 June 2021

Abstract

Background Toe systolic pressures (TSP) are an essential part of a lower limb wound assessment. They assist wound clinicians ascribe aetiology, and predict wound healing capacity, assisting with triage and ongoing care planning.

TSP method The toe cuff should be placed around the base of the hallux securely. The photoplethysmography (PPG) probe should be placed on the distal pulp of the hallux so that the probe is flush with the skin surface, affixed with tape. Wait for a strong cyclical signal from the probe to appear on the Doppler screen. Next, inflate the cuff, taking care not to overinflate. Once a flat line with the PPG signal is obtained on the screen, inflate a further 20mmHg. Then, slowly deflate the cuff at a rate of 2–5mmHg per second. The first regular cyclical waveform on the screen is the point at which the pressure in the artery is greater than the pressure in the cuff, and is considered the TSP.

Discussion TSP are a valuable measure to include in wound practice. TSP are an important part of wound staging systems, can predict wound healing capacity and help guide ongoing management.

Background

Toe systolic pressures (TSP) are an essential part of a thorough lower limb wound assessment¹. They assist wound clinicians ascribe aetiology, and predict wound healing capacity, both of which assist with triage and ongoing care planning¹. TSP are normally performed in combination with other vascular testing methods, such as pulse palpation, continuous wave Doppler waveform analysis, and ankle-brachial indices (ABI). TSP were first described by Lezack and Carter in 1970² and are now included in international guidelines for wound management³, lower extremity vascular assessment⁴, and wound grading and prediction systems such as Wifl (wound, ischaemia, foot infection)⁵.

TSP are considered advantageous compared to ABI in specific patient cohorts because digital arteries are less likely impacted by medial arterial calcification (MAC)^{6,7}. MAC is the formation of hydroxyapatite in the tunica media of the artery which leads to stiffening and loss of elasticity in the vessel and can render them incompressible⁸. This most commonly occurs in the calf vessels and is known to interfere with the accuracy of lower limb vascular testing methods, and also makes some surgical interventions more challenging⁹. In subgroups of people with diabetes, such as those with severe renal disease and advanced age (>75 years), MAC is more prevalent¹⁰⁻¹²; therefore, TSP are generally preferred as predictors of wound healing over ankle pressures in these populations¹.

TSP are generally considered to be approximately 6–10mmHg less than a brachial pressure due to the increased distance from the heart compared to brachial pressures and the smaller size of the digital arteries¹³. TSP are also used in calculations of toe-brachial indices (TBIs), which have been demonstrated to be a valuable adjunctive vascular testing method for identifying peripheral artery disease (PAD)¹⁴.

Peta Tehan* PhD, FFPM RCPS (Glasg)
Senior Lecturer Podiatry, School of Health Sciences, College of Health, Medicine and Wellbeing, BE130 Health Precinct, University of Newcastle, Ourimbah, NSW 2258 Australia
Email Peta.Tehan@newcastle.edu.au

Martin Fox FFPM RCPS (Glasg), FCPodM
Vascular Specialist Podiatrist, Manchester Leg Circulation Service, Manchester Local Care Organisation, Harpurhey Health Centre
United Kingdom

Joseph L Mills MD, FACS
Professor and Chief, Division of Vascular Surgery and Endovascular Therapy, Michael E DeBakey Department of Surgery, Baylor College of Medicine
Houston, Texas, United States of America

* Corresponding author

TSP of less than 50mmHg are generally indicative of reduced healing capacity in foot wounds, with lower values associated with worse outcomes, including higher likelihood of non-healing and amputation¹⁵. Chronic limb threatening ischaemia (CLTI) is considered when a patient presents with tissue loss, ischaemic rest pain, and a TSP of less than 60mmHg¹⁶. However, there is little data available on what a normal TSP should be, with estimates ranging from 70–100mmHg^{17,18}, and since an upper limit has not yet been established, it is currently unclear what the TSP of a digit affected by MAC would be. Generally, TSP values of greater than 140mmHg in the absence of hypertension could be considered elevated and clinicians could hold a high suspicion of MAC in the digital arteries. However, it should be noted that there is a clear lack of evidence to guide practice on this particular aspect of TSP interpretation, so interpretation of high values should be done with caution. TSP should always be correlated with the pedal and digital arterial waveforms, and interpreted in the context of other vascular testing methods, the patient, and their wound or other presenting clinical picture.

Measurement of TSP has been demonstrated to have good levels of inter- and intra-tester reliability^{19,20}. TSP are also good predictors of wound healing capacity in diabetic foot ulceration²¹ and chronic foot ulceration¹⁹, and predicting outcomes post minor amputation²². As the evidence for their use continues to grow, TSP are increasingly being recognised for their clinical utility; however, many clinicians remain uncertain on the measurement technique^{23,24}. This narrative review aims to clarify the technical aspects of TSP measurement for clinicians who manage lower limb wounds.

TSP measurement technique

Test conditions

There are a number of environmental and physical factors which need to be observed in order to optimise accuracy of TSP measurement. The room temperature should be controlled, and ideally be between 22–24°Celsius²⁵. Skin blood flow is influenced by temperature, with cold temperatures inducing vasoconstriction which will unduly influence the TSP, leading to an inaccurate reading. The patient should also be rested in a supine position for a minimum of 10 minutes prior to the measurement being taken to stabilise blood pressure²⁶. Maintaining this supine position during measurement is also important to ensure the toes are at the same level as the heart.

Other factors known to influence accuracy of measurement include smoking, caffeine intake, and exercise in the 2 hours prior to measurement²⁷, therefore these should also be avoided to optimise the accuracy of the measurement. Limb tremor, sudden movement and hyperaemia are also known to interfere with the photoplethysmography (PPG) signal²⁸. Finally, vaso-neural disorders such as Raynaud's phenomenon are known to influence TSP measurement, with structural and functional changes in the arteries²⁹ leading to potential reductions in TSP readings.

Equipment

Specialised equipment is required to measure TSP and has been cited as one of the barriers to measurement by clinicians^{23,24}. Several different types of equipment can be used to measure TSP, including PPG mercury strain gauge, laser Doppler, and continuous wave Doppler. The most commonly used is the PPG which is an infra-red probe which uses opacity changes, and changes in blood volume, to produce a waveform³⁰.

A sphygmomanometer and digital pneumatic cuff are also required, and care should be taken to choose the right sized cuff for the digit being measured. The cuff should be approximately 1.5 times the diameter of the digit. Cuffs that are too small will result in falsely high pressures and cuffs that are too large will result in erroneously low pressures³¹. Readily available cuffs come in three sizes, 1.5, 1.9 and 2.5cm, with 2.5cm suitable for most halluces and 1.5cm suitable for most other digits³². A Doppler unit that is compatible with the probe is required, and a range of equipment is currently available. Manual measures are currently preferred³³; however, some automated devices have been validated and shown to have adequate levels of reliability and accuracy³⁴. With practice, a manual TSP measurement should take less than 5 minutes to perform including the application and removal of testing equipment.

Technique

Normal hand hygiene should be followed and wearing of gloves is recommended, particularly in wounded patients. The cuff should be placed around the base of the hallux carefully, ensuring that no excess pressure is exerted on the digit whilst not inflated, however sufficiently firm to keep it in place. Care should also be taken with placement of the PPG lead to minimise any effect of lead movement or sway on the signal. In the absence of a hallux, or if a wound prevents placement of the probe, the second digit can be used as a substitute and has been shown to be equivalent³².

There is currently a lack of evidence for the accuracy of other lesser digits; however, smaller digits make cuff position and



Figure 1. PPG probe placement and PPG waveform

PPG placement challenging. In these cases, other methods of non-invasive vascular assessment can be considered, such as ankle pressures or continuous wave Doppler, or further information from vascular imaging and/or a vascular consultant could be sought, particularly if the limb appears threatened with the overall clinical presentation.

The PPG probe should be placed on the distal pulp of the hallux so that the probe is flush with the skin surface to ensure that no light can enter the underside of the probe (Figure 1). The PPG can be taped into place with hypoallergenic tape. It is important to ensure that the leads are positioned well, so as not to be moved during inflation of the cuff. The PPG probe is very sensitive to movement, so movement from the patient and also from an inflating cuff can result in interference to the waveform on the screen.

Wait for a strong cyclical signal from the probe to appear on the Doppler unit screen (Figure 1). The PPG waveform can be interpreted qualitatively, with dampening, delay and diminishment of the PPG waveform occurring with increasing severity of PAD³⁵. Next, inflate the cuff, taking care not to overinflate as the smaller sized cuff does not take much

effort to inflate and can easily burst if overinflated. Do not inflate the cuff greater than 250mmHg to avoid damaging the cuff. The waveform should generally flatten between 150–200mmHg; this may be lower in a severely diseased patient. Once a flat line with the PPG signal is obtained on the screen, inflate a further 20mmHg. Then, slowly deflate the cuff at a rate of 2–5mmHg per second. The first regular cyclical waveform (blip) on the screen (Figure 2) is the point at which the pressure in the artery is greater than the pressure in the cuff, and is considered the TSP.

This is often quite small and subtle, and commonly missed by novices. Another common error is waiting for the waveform to return to its baseline size and shape – this will result in an erroneous measure (value lower than the actual TSP). If an error is made, and the pressure needs to be re-measured, the clinician should wait 3 minutes before taking a subsequent measure. This will avoid a hyperaemic response which will artificially increase the TSP value temporarily³⁶. Once satisfied with the TSP measurement, release the sphygmomanometer pressure and remove equipment, ensuring that all equipment is decontaminated with an alcohol or neutral wipe at the end of the testing procedure. A brief step-by-step guide for TSP




Figure 2. Two examples of the return of the PPG waveform – a TSP of 105mmHg and a TSP of 19mmHg

Table 1. Brief instructions for TSP measurement

Equipment needed	Method (step-by-step)
<ul style="list-style-type: none"> • Handheld Doppler/PPG unit • Manual sphygmomanometer • Toe cuffs (2.5, 1.9 and 1.5cm) • Photoplethysmography probe (PPG) • Cleansing wipes • Hand washing facilities, hand gel and examination gloves 	<ul style="list-style-type: none"> • Gain patient consent • Perform hand hygiene • Rest the patient in a supine (flat) position for a minimum of 10 minutes • Place appropriately sized toe cuff around hallux • Attach the toe cuff to the sphygmomanometer • Attach the PPG sensor to the pulp of the hallux, affixing with hypoallergenic tape • Insert the PPG probe into the Doppler/PPG unit • Switch the unit on and wait for a strong, regular signal • Evaluate the waveform qualitatively • Inflate the toe cuff with the sphygmomanometer lightly and slowly until the waveform flattens (normally between 150–200mmHg), then inflate a further 20mmHg • Slowly release the cuff at a rate of 2–5mmHg per second • The first regular cyclical waveform is considered the toe systolic pressure • Remove equipment, decontaminate equipment, and perform hand hygiene

measurement, including equipment required, is located in Table 1.

Use of TSP in international guidelines

TSP are a useful clinical test and as such are included in multiple international guidelines for wound and vascular assessment. The International Working Group for the Diabetic Foot (IWGDF) guidelines suggest clinical vascular examination of all patients with diabetes and foot ulceration be conducted in order to identify PAD³, with TSP of less than or equal to 30mmHg necessitating urgent imaging and revascularisation. Other international guidelines, such as the American Heart Association guidelines for lower extremity PAD⁴, stipulate the use of TSP as a useful adjunct test to identify significant arterial occlusive disease. Wound grading and prediction systems such as Wifl⁵ use and prefer TSP measurement to fulfil the system's risk stratifications including risk of amputation and benefit of revascularisation. The global guidelines on CLTI mandate assessment of foot perfusion in patients with diabetic foot ulcer and/or those suspected of having CLTI (presence of PAD in combination with rest pain, gangrene or lower limb ulceration >2 weeks' duration)¹⁶. These guidelines suggest that assessment should include TSP and TBI, with TBI values of less than 0.7 considered abnormal, and TSP values less than 30mmHg typically associated with advanced ischaemia. This recommendation was made to ensure that CLTI is not missed in such patients for whom a foot ulcer may be a precursor to amputation. As the evidence for TSP and TBI continue to grow, the uptake into clinical guidelines will also likely increase.

Care planning and onward referral

Ideally, foot wound and vascular care would be provided as part of an integrated, multidisciplinary limb salvage team, but this is not possible or practical in all clinical settings. This is particularly the case in areas of rural and remote Australia, where specialised integrated services are lacking. Pulse examination alone is necessary but insufficient to evaluate perfusion in the wounded lower limb⁹. A validated limb-staging system such as Wifl should be used to stage the limb and its predicted risk of amputation, as well as to guide initial therapy. In patients with active foot wounds, the preferred measurement of foot perfusion is TSP, although ABIs may be used in situations where they can be expected to be reliable, such as in patients who have compressible vessels without MAC.

Wifl staging is based on grading the wound, the perfusion, and the presence and severity of wound infection. The combination of grades in each of these three domains is used to place the patient into clinical stages (1 to 4) that correlate with 1-year amputation risk. Details of this system are available in the original publications^{5,16} and a free mobile application is available to calculate Wifl grades and clinical stages. Any patient with a foot ulcer and ischaemia Grade 2

or 3, as well as any patient with Wifl Clinical Stage ≥ 2 should be referred urgently to a vascular surgeon/specialist for further assessment of PAD. Where the wound assessment takes place in remote areas, locations far away from vascular departments, during times such as the recent COVID pandemic, or when the clinician is unsure of severity, a timely telemedicine consultation with the vascular team, including any Wifl results on hand, may be helpful to assist decision-making, triage, and avoidance of amputations.

Clinical Stage 1 patients, and those patients with Grade 0 ischaemia based on TSP >60mmHg, will usually heal with offloading and optimal wound care. However, even in these cases, vascular referral should be considered if wound healing stalls, or fails to occur, after 4–6 weeks of optimal care. In addition, it should be noted that TSP reflects toe and forefoot perfusion. There are patients with regional malperfusion of the foot, so these measurements should be interpreted with caution in patients with midfoot and hindfoot ulcers, particularly if the posterior tibial artery may be severely stenosed or occluded³⁷, which may present as a low posterior tibial artery systolic pressure or an absent or a dampened monophasic Doppler waveform signal³⁸.

In people with chronic lower limb oedema, ABI has long been suggested as part of the holistic assessment to determine safety prior to considering the use of compression bandaging or hosiery³⁹, with most guidelines stipulating an ABI of greater than or equal to 0.8 is required to proceed with high-level compression⁴⁰. TBI has also been advocated where ABI is difficult due to the presence of oedema or in cases of MAC⁴¹. Despite these recommendations, TBI seems to be rarely used by lower limb clinicians to date, despite barriers to ABI performance, such as limb oedema, being common in people admitted in hospital for foot ulceration⁴². Despite little published research to date in this area, TSP may therefore have an important role in the assessment and management of chronic limb oedema, particularly where ABI is unsuitable for the reasons discussed above.

Conclusion

TSP are a valuable vascular testing method to include in wound assessment practice, particularly in patients with active foot ulceration, diabetes or renal disease. TSP are an important part of wound staging systems, and can predict wound healing capacity and help guide ongoing management. TSP play an important part of a comprehensive wound and lower limb vascular assessment and, once learned, are quick and easy to perform.

Conflicts of interest

None of the authors have any conflicts of interest to declare.

Ethics statement

Not applicable.

Funding

The authors received no funding for this study.

References

1. Tehan PE, Barwick AL, Casey SL, Lanting SM, Chuter VH. Accurate noninvasive arterial assessment of the wounded lower limb: a clinical challenge for wound practitioners. *Int J Lower Extrem Wounds* 2020;19(3):215–226.
2. Lezack J, Carter S. Systolic pressures in the extremities of man with special reference to the toes. *Canadian J Physiol Pharmacol* 1970;48(7):469–474.
3. Hinchliffe RJ, Forsythe RO, Apelqvist J, Boyko EJ, Fitrige R, Hong JP, Katsanos K, Mills JL, Nikol S, Reekers J et al. Guidelines on diagnosis, prognosis, and management of peripheral artery disease in patients with foot ulcers and diabetes (IWGDF 2019 update). *Diabetes/Metabolism Res Rev* 2020;36 Suppl 1:e3276.
4. Rooke TW, Hirsch AT, Misra S, Sidawy AN, Beckman JA, Findeiss LK, Golzarian J, Gornik HL, Halperin JL, Jaff MR et al. 2011 ACCF/AHA focused update of the guideline for the management of patients with peripheral artery disease (updating the 2005 guideline): a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. *J Am College Cardiol* 2011;58(19):2020–2045.
5. Mills JL, Sr., Conte MS, Armstrong DG, Pomposelli FB, Schanzer A, Sidawy AN, Andros G. The Society for Vascular Surgery lower extremity threatened limb classification system: risk stratification based on wound, ischemia, and foot infection (Wlfi). *J Vasc Surg* 2014;59(1):220–234.e222.
6. Leskinen Y, Salenius JP, Lehtimäki T, Huhtala H, Saha H. The prevalence of peripheral arterial disease and medial arterial calcification in patients with chronic renal failure: requirements for diagnostics. *Am J Kidney Dis* 2002;40(3):472–479.
7. Young M, Adams JE, Anderson GF, Boulton AJM, Cavanagh PR. Medial arterial calcification in the feet of diabetic patients and matched non-diabetic control subjects. *Diabetologia* 1993;36:615–621.
8. Ho CY, Shanahan CM. Medial arterial calcification: an overlooked player in peripheral arterial disease. *Arterioscler, Thromb Vasc Biol* 2016;36(8):1475–1482.
9. Hinchliffe RJ, Brownrigg JR, Apelqvist J, Boyko EJ, Fitrige R, Mills JL, Reekers J, Shearman CP, Zierler RE, Schaper NC. IWGDF guidance on the diagnosis, prognosis and management of peripheral artery disease in patients with foot ulcers in diabetes. *Diabetes Metab Res Rev* 2016;32 Suppl 1:37–44.
10. Suominen V, Rantanen T, Venermo M, Saarinen J, Salenius J. Prevalence and risk factors of PAD among patients with elevated ABI. *Eur J Vasc Endovasc Surg* 2008;35(6):709–714.
11. Potier L, Halbron M, Bouilloud M, Dadon M, Le Doeuff J. Ankle-to-brachial ratio index underestimates the prevalence of peripheral occlusive disease in diabetic patients at high risk for arterial disease. *Diabetes Care* 2009;32(4):e44.
12. Emanuele M, Buchanan B, Abaira C. Elevated leg systolic pressures and arterial calcification in diabetic occlusive vascular disease. *Diabetes Care* 1981;4(Mar-Apr):289–292.
13. Nielsen P, Bell G, Lassen NA. The measurement of digital systolic blood pressure by strain gauge technique. *Scand J Clin Lab Invest* 1972;29:371.
14. Tehan PE, Bray A, Chuter VH. Non-invasive vascular assessment in the foot with diabetes: sensitivity and specificity of the ankle brachial index, toe brachial index and continuous wave Doppler for detecting peripheral arterial disease. *J Diabetes Complicat* 2016;30(1):155–160.
15. Sonter J, Ho A, Chuter VH. The predictive capacity of toe blood pressure and the toe-brachial index for foot wound healing and amputation: a systematic review and meta-analysis. *Wound Practice & Res* 2014;22(4):208–217.
16. Conte MS, Bradbury AW, Kolh P, White JV, Dick F, Fitrige R, Mills JL, Ricco J-B, Suresh KR, Murad MH. Global vascular guidelines on the management of chronic limb-threatening ischemia. *Eur J Vasc Endovasc Surg* 2019;58(1):S1-S109, e133.
17. Tehan PE, Barwick AL, Sebastian M, Chuter VH. Diagnostic accuracy of resting systolic toe pressure for diagnosis of peripheral arterial disease in people with and without diabetes: a cross-sectional retrospective case-control study. *J Foot Ankle Res* 2017;10:58.
18. Patry J, Laurencelle L, Bélisle J, Beaumier M. Vascular assessment in patients with a lower limb wound: a correlational study of photoplethysmography and laser Doppler flowmetry toe pressure techniques. *J Diabetes Sci Technol* 2020; <https://doi.org/10.1177/1932296820979973>.
19. Sonter J, Chuter V, Casey S. Intra-tester and inter-tester reliability of toe pressures and toe brachial indices: a systematic review Manuscript submitted for publication; in press.
20. Scanlon C, Park K, Mapletoft D, Begg L, Burns J. Interrater and intrarater reliability of photoplethysmography for measuring toe blood pressure and toe-brachial index in people with diabetes mellitus. *J Foot Ankle Res* 2012;5(13).
21. Forsythe RO, Apelqvist J, Boyko EJ, Fitrige R, Hong JP, Katsanos K, Mills JL, Nikol S, Reekers J, Venermo M, et al. Performance of prognostic markers in the prediction of wound healing or amputation among patients with foot ulcers in diabetes: a systematic review. *Diabetes/Metabolism Res Rev* 2020;36 Suppl 1:e3278.
22. Linton C, Searle A, Hawke F, Tehan PE, Sebastian M, Chuter V. Do toe blood pressures predict healing after minor lower limb amputation in people with diabetes? A systematic review and meta-analysis. *Diabetes Vasc Dis Res* 2020;17(3): doi: 10.1177/1479164120928868.
23. Tehan PE, Fox M, Stewart S, Matthews S, Chuter VH. Lower limb vascular assessment techniques of podiatrists in the United Kingdom: a national survey. *J Foot Ankle Res* 2019;12(1):31.
24. Tehan PE, Chuter VH. Vascular assessment techniques of podiatrists in Australia and New Zealand: a web-based survey. *J Foot Ankle Res* 2015;8(1):71.
25. Sawka A, Carter SA. Effect of temperature on digital systolic pressures in lower limb in arterial disease. *Circulation* 1992;85:1097–1101.
26. Sadler S, Hawke F, Chuter V. Effect of pre-rest test duration on toe and ankle systolic pressure measurements. *J Foot Ankle Res* 2013;6(Suppl 1):13.
27. Campbell N, Chockalingham A, Fodor JG, McKay DW. Accurate, reproducible measurement of blood pressure. *Canadian Med Assoc J* 1990;143(1):19–24.
28. Pérez-Martin A, Meyer G, Demattei C, Böge G, Laroche J-P, Quéré I, Dauzat M. Validation of a fully automatic photoplethysmographic device for toe blood pressure measurement. *Eur J Vasc Endovasc Surg* 2010;40(4):515–520.
29. Ho M, Belch JFF. Raynaud's phenomenon: state of the art 1998. *Scand J Rheumatol* 1998;27(5):319–322.
30. Rose SC. Noninvasive vascular laboratory for evaluation of peripheral arterial occlusive disease: Part I: hemodynamic principles and tools of the trade. *J Vasc Interv Radiol* 2000;11(9):1107–1114.

31. Pålsson H-I, Laskar C, Stark K, Andersson A, Jogestrand T, Wahlberg E. The optimal cuff width for measuring toe blood pressure. *Angiol* 2007;58(4):472–476.
32. Bhamidipaty V, Dean A, Yap SL, Firth J, Barron M, Allard B, Chan ST. Second toe systolic pressure measurements are valid substitutes for first toe systolic pressure measurements in diabetic patients: a prospective study. *Eur J Vasc Endovasc Surg* 2015;49(1):77–82.
33. The National Institute for Health and Care Excellence (NICE). Lower limb peripheral arterial disease: diagnosis and management (CG147); 2020. Available from: <http://www.nice.org.uk/guidance/CG147>
34. Sonter J, Tehan P, Chuter V. Toe brachial index measured by automated device compared to duplex ultrasonography for detecting peripheral arterial disease in older people. *Vascular* 2017;25(6):612–617.
35. Allen J. Photoplethysmography and its application in clinical physiological measurement. *Physiol Measure* 2007;28(R1).
36. Cracowski J-L, Minson CT, Salvat-Melis M, Halliwill JR. Methodological issues in the assessment of skin microvascular endothelial function in humans. *Trends Pharmacol Sci* 2006;27(9):503–508.
37. Montero-Baker M, Rocha B, Armstrong D, Goshima K, Miller J. Diagnosis and endovascular management of segmental heel ischemia. *Int J Clin Cardiol* 2018;5:117.
38. Kim ES, Sharma AM, Scissons R, Dawson D, Eberhardt RT, Gerhard-Herman M, Hughes JP, Knight S, Marie Kupinski A, Mahe G. Interpretation of peripheral arterial and venous Doppler waveforms: a consensus statement from the Society for Vascular Medicine and Society for Vascular Ultrasound. *Vasc Med* 2020;25(5):484–506.
39. Kelechi TJ, Johnson JJ, Yates S. Chronic venous disease and venous leg ulcers: an evidence-based update. *J Vascular Nurs* 2015;33(2):36–46.
40. Weller CD, Team V, Ivory JD, Crawford K, Gethin G. ABPI reporting and compression recommendations in global clinical practice guidelines on venous leg ulcer management: a scoping review. *Int Wound J* 2019;16(2):406–419.
41. Doherty D, Morgan P, Moffatt C. Hosiery in lower limb lymphoedema. *J Lymph* 2009;1:30–37.
42. Quéré I, Palmier S, Nøerregaard S, Pastor J, Sykorova M, Dring E, Franks PJ, Murray S, Keeley V, Bermark S. LIMPRINT: estimation of the prevalence of lymphoedema/chronic oedema in acute hospital in in-patients. *Lymphat Res Biol* 2019;17(2):135–140.